

Booster Buddy Seat Replacement Strap Instructions

P/N 289-KIT

Tools Required:

Phillips Head Screwdriver

Materials Supplied:

Replacement Strap

Two 8-16 x 1/2" Plastite Thread Rolling Screws

1. To remove the strap for replacement, unscrew the chafe tabs on each side of the booster using a Phillips head screwdriver. Discard old screws and strap.

2. To replace the new strap, ensure that the flat part of each chafe tab is against the plastic. Using the screws provided, screw through the hole in the chafe tabs and into the plastic.

